



Recombinant Human HIRA-interacting protein 3 (HIRIP3)

Product Code	CSB-YP883612HU
Abbreviation	HIRIP3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BW71
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAREKEMQEF TRSFFRGRPD LSTLTHSIVR RRYLAHSGRS HLEPEEKQAL KRLVEEELLK MQVDEAASRE DKLDLTKKGK RPPTPCSDPE RKRFRFNSES ESGSEASSPD YFGPPAKNGV AAEVSPAKEE NRRASKAVE ESSDEERQRD LPAQRGEES EEEKGYKKGK TRKKPVVKKQ APGKASVSRK QAREESESE AEPVQRTAKK VEGNKGTKSL KESEQESEEE ILAQKKEQRE EEVEEEEEKEE DEEKGDWKPR TRSNGRRKSA REERSCKQKS QAKRLLGDSD SEEEQKEAAS SGDDSGRDRE PPVQRKSEDR TQLKGGKRLS GSSEDEEDSG KGEPTAKGSR KMARLGSTSG EESDLEREVS DSEAGGGPQG ERKNRSSKKS SRKGRTRSSS SSSDGSPEAK GGKAGSRRG EDHPAVMRLK RYIRACGAHR NYKLLGSCC SHKERLSILR AELEALGMKG TPSLGKCRAL KEQREEAAEV ASLDVANIIS GSGRPRRRTA WNPLGEAAPP GELYRRTLDS DEERPRPAPP DWSHMRGIIS SDGESN
Source	Yeast
Target Names	HIRIP3
Protein Names	Recommended name: HIRA-interacting protein 3
Expression Region	1-556
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The HIRA protein shares sequence similarity with Hir1p and Hir2p, the two corepressors of histone gene transcription characterized in the yeast, <i>Saccharomyces cerevisiae</i> . The structural features of the HIRA protein suggest that it may function as part of a multiprotein complex. Recently, several cDNAs encoding HIRA-interacting proteins, or HIRIPs, have been identified. In vitro, the HIRIP3 gene product binds HIRA, as well as H2B and H3 core histones, indicating that a complex containing HIRA-HIRIP3 could function in some



aspects of chromatin and histone metabolism.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.