



Recombinant Human Testis-expressed sequence 13A protein (TEX13A)

Product Code	CSB-BP883626HU
Abbreviation	TEX13A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BXU3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MALRPEDPSS GFRHSNVVAF INEKMARHTK GPEFYLENIS LSWEKVEDKL RAILEDSEVP SEVKEACTWG SLALGVRF AH RQAQLQRHRV RWLHGFAKLH KSAAQALASD LKKLREQQET ERKEAASRLR MAQTSLVEVQ KERDKELVSP HEWEQGAGWP GLATAGGVCT EGAAEEEEEEA AVAAAGAAGG KGAEQQRDV EVVAAPVEAM APPVEAGAAP METQFPHVEA RAASMETTEK LERILLQLLG DADQEKYTYW GQKEGDLRSV ETATSYFSGT TNPWSRASSE PLPVQLPASY SYSYSSPFSS FSDIPTISPP QATVTAPVPP QLPSDWEAFD TSLWSDGGPH RIDHQEHPRD RRYSEPHQQR PPVYRRPGDW DCPWCNAVNF SRRDTCFDCG KGIWLQKPH
Source	Baculovirus
Target Names	TEX13A
Protein Names	Recommended name: Testis-expressed sequence 13A protein
Expression Region	1-409
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene is similar to a mouse gene that is expressed in the testis.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.