



# Recombinant Human Pre-B-cell leukemia transcription factor 4 (PBX4)

<b>Product Code</b>	CSB-MP883638HU
<b>Abbreviation</b>	PBX4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9BYU1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAAPPRPAPS PPAPRRLDTS DVLQQIMAIT DQSLDEAQAR KHALNCHRMK PALFSVLCEI KEKTVVSIRG IQDEDPPDAQ LLRLDNMLLA EGVC RPEKRG RGGAVARAGT ATPGGCPNDN SIEHSDYRAK LSQIRQIYHS ELEKYEQACR EFTTHVTNLL QEQRMRPVS PKEIERMVG A IHGKFSAIQM QLKQSTCEAV MTLRSRLDA RRKRRNFSKQ ATEVLNEYFY SHLNNPYPSE EAKEELARKG GLTISQVSNW FGNKRIRYKK NMGKFQEEAT IYTGKTAVD T TEVGVPGNHA SCLSTPSSGS SGPFLPSAG DAFLTLRTLA SLQPPPGGGC LQSQAQGSWQ GATPQPATAS PAGDPGSINS STSN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PBX4
<b>Protein Names</b>	Recommended name: Pre-B-cell leukemia transcription factor 4 Alternative name(s): Homeobox protein PBX4
<b>Expression Region</b>	1-374
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a homeodomain protein with similarity to a transcription factor involved in translocations in pre-B-cell leukemias. The presence of a homeobox domain suggests that this protein acts as a transcription factor, but this function has not been demonstrated.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.