



Recombinant Human Rho-related BTB domain-containing protein 2 (RHOBTB2)

Product Code	CSB-YP883640HU
Abbreviation	RHOBTB2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BYZ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSMDYERP NVETIKCVVV GDNVAVGKTRL ICARACNATL TQYQLLATHV PTVWAIQYR VCQEVLESLR DVVDDVSVSL RLWDTFGDHH KDRRFAYGRS DVVLCFSIA NPNSLHHVKT MWYPEIKHFC PRAPVILVGC QLDLRYADLE AVNRARRPLA RPIKPNEILP PEKGREVAKE LGIPYYETSV VAQFGIKDVF DNAIRAALIS RRHLQFWKSH LRNVQRPLLQ APFLPPKPPP PIIVVPDPPS SSECPAHL EDPLCADVIL VLQERVIFA HKIYLSTSSS KFYDLFLMDL SEGELGGPSE PGGTHPEDHQ GHSDQHSHHH HHHHGRDFLL RAASFDVCES VDEAGGSGPA GLRASTSDGI LRGNGTGYLP GRGRVLSSWS RAFVSIQEEM AEDPLTYKSR LMVVVKMDSS IQPGPFRAVL KYLYTGELDE NERDLMHIAH IAELLEVFDL RMMVANILNN EAFMNQEITK AFHVRRTNRV KECLAKGTFS DVTFILDDGT ISAHKPLLIS SCDWMAAMFG GPFVESSTRE VVPYTSKSC MRVLEYLYT GMFTSSPDLD DMKLIILANR LCLPHLVALT EQYTVTGLME ATQMMVDIDG DVLVFLELAQ FHCAYQLADW CLHHICTNYN NVCCKFPRDM KAMSPENQEY FEKHRWPPVW YLKEEDHYQR ARKEREKEDY LHLKRQPKRR WLFWNSSPSS SSSAASSSSP SSSSAVV
Source	Yeast
Target Names	RHOBTB2
Protein Names	Recommended name: Rho-related BTB domain-containing protein 2 Alternative name(s): Deleted in breast cancer 2 gene protein p83
Expression Region	1-727
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.