



Recombinant Human Uncharacterized protein KIAA1751 (KIAA1751)

Product Code	CSB-EP883654HU-B
Abbreviation	KIAA1751
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C0B2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEDDGSELLPE DELLADALLL EDERDELEDP EFDIKCLLQE AEDDVDPGHS SSVKELDTDA DKLKKKTAED RTQAFHLRQN LSALDKMHEE QELFTEKMRG ELRACRQRD LIDKQQEAVA AEIATEEEAG NMAAVGRLQA VSRRLFAELE NERDLQSRTE AVLKESENTM WHIEIQEGRLEAFRTADREE VEATGRRLQV RAAEQLCREQ EALGKVERNR LLRIRKSLNT QKELGLRHQK LLEDARKNHK VAVRFLKASL GRIREQEKKE EMECHEYMRR RMDAVVALKG SISRNRDTRL KFQAWDRAKA ELAEQRVQAE KKAILAQGRD AFRHLVHQRR RQELEAQKRA FEEEQKLRKQ EIISRILKEE AEEEKRKKQH PPTSARHRLT LRDKTWNYIS DFCKKTTVPT NTYTLDYEAA AGPGPSRLLE VVSSELIQGD PGASSEETL AEPEISGLWN EDYKPYQVPK EDVDRKPVGG TKMDKDILER TVERLRSRVV HKQVWVWVREF QGRPFNSKPE LLHFQDFDIG KVKYKKITLV NTTYTINYCK LVGVVEHLRD FIHVDFDPPG PLSAGMSCEV LVTFKPMINK DLEGNISFLA QTGEFSVPLK CSTKKCSLSL DKELIDFGSY VVGETTSRTI TLTNVGGLGT TFKFLPASEP CEMDDSQSAL KLSLLTYED KSLYDKAATS FSEQQLEGTE SSQADMQSRLK ELEKLDKEQE EEQPAGEPGC QAKPGWAGGE GLRAAFQGAE GPRGGLESRE FS
Source	E.coli
Target Names	CFAP74
Protein Names	Recommended name: Uncharacterized protein KIAA1751
Expression Region	1-762
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.