



Recombinant Human Kelch-like protein 4 (KLHL4)

Product Code	CSB-MP883657HU
Abbreviation	KLHL4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C0H6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSVSGKKEFD VKQILRLRWR WFSHPFQGST NTGSCLQQEG YEHRGTPVQG RLKSHSRDRN GLKKSNSPVH HNILAPVPGP APAHQRAVQN LQQHNLIHF QANEDTPKSV PEKNLFKEAC EKRAQDLEMM ADDNIEDSTA RLDTQHSEDM NATRSEEQFH VINHAEQTLR KMENYLKEKQ LCDVLLIAGH LRIPAHLVL SAVSDYFAAM FTNDVLEAKQ EEVRMEGVDP NALNSLVQYA YTGVLQLKED TIESLLAAC LLQLTQVIDV CSNFLIKQLH PSNCLGIRSF GDAQGCTELL NVAHKYTMEH FIEVIKNQEF LLLPANEISK LLCSDDINVP DEETIFHALM QWVGHDVQNR QGELGMLLSY IRLPLLPPQL LADLETSSMF TGDLECQKLL MEAMKYHLLP ERRSMMQSPR TKPRKSTVGA LYAVGGMDAM KGTTTIEKYD LRTNSWLHIG TMNGRRLQFG VAVIDNKLYV VGGRDGLKTL NTVECFNPVG KIWTVMP PMS THRHGLGVAT LEGPMYAVGG HDGWSYLNTV ERWDPEGRQW NYVASMSTPR STVGVALNN KLYAIGGRDG SSCLKSMEYF DPHTNKWSLC APMSKRGGV GVATYNGFLY VVGHDAPAS NHCSRLSDCV ERYDPKGDSW STVAPLSVPR DAVAVCPGLD KLYVVGGYDG HTYLNTVESY DAQRNEWKEE VPVNIGRAGA CVVVVKLP
Source	Mammalian cell
Target Names	KLHL4
Protein Names	Recommended name: Kelch-like protein 4
Expression Region	1-718
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a member of the kelch family of proteins, which are characterized by kelch repeat motifs and a POZ/BTB protein-binding domain. It is thought that kelch repeats are actin binding domains. However, the specific function of this protein has not been determined. Alternative splicing of this gene results in two transcript variants encoding different isoforms.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.