



# Recombinant Glyoxylate reductase (gyaR)

<b>Product Code</b>	CSB-EP883676TMQ
<b>Abbreviation</b>	gyaR
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9C4M5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Thermococcus litoralis (strain ATCC 51850 / DSM 5473 / JCM 8560 / NS-C)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKPKVFITRQ IPENGIKMIE KFYEIELWKD PKAPPRGVLL EKVREVDALV TLVTDKVDKE LLENAPKLLKI IAQYAVGYDN IDIEEATKRG IYVTNTPGVL TDATADLafa LLLAVARRIV EADAFVRSGE WKKSEVGWHP LMFLGYGLKG KTLGIVGFGR IGQALAKRAK GFGMKIIYYS RTRKPEAEeee IGAEYVDFET LLKESDFISL HVPLTKETYH MIGEkelKLM KPNAILINTS RGAVVDTNAL IKALKEGWIA GAGLDVFEEe PYyNEELFKL KNVVLAPHIG SATHEAREGM AELVAKNLIA FAKGEIPPNL VNKDVLtSSP P
<b>Source</b>	E.coli
<b>Target Names</b>	gyaR
<b>Protein Names</b>	Recommended name: Glyoxylate reductase EC= 1.1.1.26
<b>Expression Region</b>	1-331
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.