



Recombinant Arabidopsis thaliana Putative F-box/FBD/LRR-repeat protein At1g66290 (At1g66290)

Product Code	CSB-EP883713DOA-B
Abbreviation	At1g66290
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9C8Y8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MDEDGERRVR TKRSCSPESS DNGSGDEV DW ISDLPEALIV LVLLNLPTKD VIKTSVLSTK WRNIWRYVPR LDLDNRH FTE FDSLASFIDR FMRSNRKVNK FKLRCGSDLD GDVDLATCSW KWIYMAIKRK VQHIDVTWPG VRIYPEIYNC ESLVSLKLSE VTLTKPEFVS LPSLKVLD WVEFYNEFAF DMLMSGCPVL ESFTLCRRYL DNVRVLRVIS QSLLSFNYG SSSKSFRRDDL VVILDAPKLE KLKLSHQLTA SFIIENLSSL VEADIDIEFN FCRGKKFDPD DLPKREMIRN FLVGLSRVKT MTIAACTLEV IYDYSRCEPL PLFPNLSFFS VEFYQKRWEI LPFFFKSCSN LKSLVLESDY FPRKRTSIIIS EPRCLLSSLE YVKIEFALDK GKMELVRYLL ENSPILKKLT LSLDHSSRKK SCVILRELIT IPRCSTSCRV IVL
Source	E.coli
Target Names	At1g66290
Protein Names	Recommended name: Putative F-box/FBD/LRR-repeat protein At1g66290
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.