



# Recombinant Arabidopsis thaliana UDP-glycosyltransferase 89B1 (UGT89B1)

<b>Product Code</b>	CSB-MP883716DOA
<b>Abbreviation</b>	UGT89B1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9C9B0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKVNEENKPK TKTHVLIFPF PAQGHMIPLL DFTHRLALRG GAALKITVLV TPKNLPFLSP LLSAVVNI EP LILPFP SHPS IPSGVENVQD LPPSGFPLMI HALGNLHAPL ISWITSHPS PVAIVSDFFL GWTKNLGIPR FDFSPSAAIT CCILNTLWIE MPTKINEDDD NEILHF PKIP NCPKYRFDQI SLYRSYVHG DPAWEFIRDS FRDNV ASWGL VVNSFTAMEG VYLEHLKREM GHDRVWAVGP IIPLSGDNRG GPTS SVSDHV MSWLDAREDN HVVYVCFGSQ VVLTKEQTLA LASGLEKSGV HFIWAVKEPV EKDSTRGNIL DGFDDR VAGR GLVIRGWAPQ VAVLRHRAVG AFLTHCGWNS VVEAVVAGVL MLTWPMRADQ YTDASLVVDE LKVGVRACEG PDTVPDPDEL ARVFADSVTG NQTERIKAVE LRKAALDAIQ ERGSSVNDLD GFIQHVVSLG LNK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	UGT89B1
<b>Protein Names</b>	Recommended name: UDP-glycosyltransferase 89B1 EC= 2.4.1.- Alternative name(s): Flavonol 3-O-glucosyltransferase UGT89B1 EC= 2.4.1.91 Flavonol 7-O-glucosyltransferase UGT89B1
<b>Expression Region</b>	1-473
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.