



# Recombinant Lactococcus lactis subsp. lactis Enolase 1 (eno1)

<b>Product Code</b>	CSB-YP883798LNG
<b>Abbreviation</b>	eno1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9CHS7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Lactococcus lactis subsp. lactis (strain IL1403) (Streptococcus lactis)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSIITDIYAR EVLDSRGNPT LEVEVYTEDG AFGRGMVPSG ASTGEHEAVE LRDGDKSRYN GLGTQKAVDN VNNIIAEAI GYEVDQQA DRAMIALDGT ENKGKLGANA ILGVSIAAAR AADELGVPL YNYLGGFNAK VLPTPMMNII NGGSHSDAPI AFQEFMIVPV GAPTFKEALR WGAEIFHALK KILKARGLET AVGDEGGFAP KFDGTEGVE TILKAIEAAG YKAGEDGVM GFDCASSEFY ENGVYDYTKF EGEKKLSA SEQVDYLEEL VSKYPIITIE DGMDENDWDG WKILTERLGK KVQLVGDDFF VTNTKYLARG IRENASNAIL IKVNQIGTLT ETFEAIEMAK EAGFTAIVSH RSGETEDSTI SDIAVATNAG QIKTGSLSRT DRMAKYNQLL RIEDQLAEVA QYKGLKAFYN LKK
<b>Source</b>	Yeast
<b>Target Names</b>	eno1
<b>Protein Names</b>	Recommended name: Enolase 1 EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase 1 2-phosphoglycerate dehydratase 1
<b>Expression Region</b>	1-433
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.