



# Recombinant Mouse Neurotensin/neuromedin N (Nts)

<b>Product Code</b>	CSB-EP883989MO-B
<b>Abbreviation</b>	Nts
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9D3P9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	SDSEEDVR ALEADLLTNM HTSKISKASP PSWKMTLLNV CSLINNVNSP AEEAGDMHDD DLVGKRKLPL VLDGFSLEAM LTIFQLQKIC RSRAFQHWI IQEDILDNVN DKNEKEEVK RKIPYIL
<b>Source</b>	E.coli
<b>Target Names</b>	Nts
<b>Protein Names</b>	Recommended name: Neurotensin/neuromedin N Cleaved into the following 4 chains: 1. Large neuromedin N Alternative name(s): NmN-125 Neuromedin N Short name= NN Short name= NmN Neurotensin Short name= NT Tail pep
<b>Expression Region</b>	23-147
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a common precursor for two peptides, neuromedin N and neurotensin. Neurotensin is a secreted tridecapeptide, which is widely distributed throughout the central nervous system, and may function as a neurotransmitter or a neuromodulator. It may be involved in dopamine-associated pathophysiological events, in the maintenance of gut structure and function, and in the regulation of fat metabolism. Tissue-specific processing may lead to the formation in some tissues of larger forms of neuromedin N and neurotensin. The large forms may represent more stable peptides that are also biologically active.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.