



# Recombinant Mouse NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial (Ndufs7)

<b>Product Code</b>	CSB-EP884070MO
<b>Abbreviation</b>	Ndufs7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9DC70
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	PSPSP SPSLSSTQSA VSKAGAGAVV PKLSHLPRSR AEYVVTKLDD LINWARRSSL WPMTFGLACC AVEMMHMAAP RYDMDFRGVV FRASPRQADV MIVAGTLTNK MAPALRKVYD QMPEPRYVVS MGSCANGGGY YHYSYSVVRG CDRIVPVDIY VPGCPPTAEA LLYGILQLQR KIKREQKLKI WYRR
<b>Source</b>	E.coli
<b>Target Names</b>	Ndufs7
<b>Protein Names</b>	Recommended name: NADH dehydrogenase [ubiquinone] iron-sulfur protein 7, mitochondrial EC= 1.6.5.3 EC= 1.6.99.3 Alternative name(s): Complex I-20kD Short name= CI-20kD NADH-ubiquinone oxidoreductase 20 kDa subunit
<b>Expression Region</b>	36-224
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a protein that is a subunit of one of the complexes that forms the mitochondrial respiratory chain. This protein is one of over 40 subunits found in complex I, the nicotinamide adenine dinucleotide (NADH):ubiquinone oxidoreductase. This complex functions in the transfer of electrons from NADH to the respiratory chain, and ubiquinone is believed to be the immediate electron acceptor for the enzyme. Mutations in this gene cause Leigh syndrome due to mitochondrial complex I deficiency, a severe neurological disorder that results in bilaterally symmetrical necrotic lesions in subcortical brain regions.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.