



Recombinant Mouse F-box/LRR-repeat protein 12 (Fbxl12)

Product Code	CSB-EP884125MO
Abbreviation	Fbxl12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9EPX5
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MATLFDLPDL VLLIEFSYLP VRDRIRISRV CHRWKRLVDD RWLWRHVDLT LYTMRPKVMW HLLRRYMASR LYSLRMGGYL FSGSQAPQLS PALMRALGQK CPNLKRLCLH VADLSMVPIT SLPSTLRTLE LHSCEISMIW LQKEQDPTVL PLLECIVLDR VPAFRDEHLQ GLTRFRALRS LVLGGTYRVT ETGLDASLQE LSYLQRLEVL GCTLSADSTL LAISRHLRDV RKIRLTVGGL SAQGLVFLEG MPVLESLECFQ GPLITPDMPT PTQIVSSCLT MPKLRVLEVQ GLGWEGQEAE KILCKGLPHC IVIVRACPKE SMDWWM
Source	E.coli
Target Names	Fbxl12
Protein Names	Recommended name: F-box/LRR-repeat protein 12 Alternative name(s): F-box and leucine-rich repeat protein 12 F-box protein FBL12
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.