



# Recombinant Rat C-terminal-binding protein 2 (Ctbp2)

<b>Product Code</b>	CSB-EP884130RA-B
<b>Abbreviation</b>	Ctbp2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9EQH5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MALVDKHKVK RQRLDRICEG IRPQIMNGPL HPRPLVALLD GRDCTVEMPI LKDLATVAFC DAQSTQEIHE KVLNEAVGAM MYHTITLTRE DLEKFKALRV IVRIGSGYDN VDIKAAGELG IAVCNIPSAA VEETADSTVC HILNLYRRNT WLYQALREGT RVQSVEQIRE VASGAARIRG ETLGLIGFGR TGQAVAVRAK AFGFSVIFYD PYLQDGIERS LGVQRVYTLQ DLLYQSDCVS LHCNLNEHNN HLINDFTIKQ MRQGAFLVNA ARGGLVDEKA LAQALKEGRI RGAALDVHES EPFSFAQGPL KDAPNLICTP HTAWYSEQAS LEMREAAATE IRRAITGRIP ESLRNCVNKE FFVTSTPWSV IDQQAIHPEL NGATYRYPPG IVGVAPGGPL PAMEGIIPGG IPVTHNLPTV AHPSQAPSPN QPTKHGDNRE HPNEQ
<b>Source</b>	E.coli
<b>Target Names</b>	Ctbp2
<b>Protein Names</b>	Recommended name: C-terminal-binding protein 2 Short name= CtBP2
<b>Expression Region</b>	1-445
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene produces alternative transcripts encoding two distinct proteins. One protein is a transcriptional repressor, while the other isoform is a major component of specialized synapses known as synaptic ribbons. Both proteins contain a NAD <sup>+</sup> binding domain similar to NAD <sup>+</sup> -dependent 2-hydroxyacid dehydrogenases. A portion of the 3 untranslated region was used to map this gene to chromosome 21q21.3; however, it was noted that similar loci elsewhere in the genome are likely. Blast analysis shows that this gene is present on chromosome 10.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.