



Recombinant Mouse Beta-parvin (Parvb)

Product Code	CSB-EP884150MO
Abbreviation	Parvb
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q9ES46
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<pre>MSSAPPRSPT PRAPKMKKDE SFLGKLGGLT ARKKKTREVT DLQEEGKSAI NSPMAPALVD IHPEDTQLEE NEERTMIDPT SREDPKFKEL VKVLLDWIND VLAERIIVK QLEEDLYDGQ VLQKLEKLA HCKLNVAEVT QSEIGQKQKL QTVLEAVQDL LRPHGWPLRW NVDSIHGKNL VAILHLLVSL AMHFRAPIHL PEHVTVQVVV VRKREGLLHS SHISEELTTT TEIMMGRFER DAFDTLFDHA PDKLNLVKKS LITFVNKHLN KLNLEVTDL E TQFADGVYLV LLLGLLEDYF VPLHNFYLT P DSFDQKVHNV AFAPFELMLDG GLKKPKARPE DVVNLDLKST LRVLYTLFTK YKDVE</pre>
Source	E.coli
Target Names	Parvb
Protein Names	Recommended name: Beta-parvin
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>