



# Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At5g42450, mitochondrial (PCMP-E102)

<b>Product Code</b>	CSB-EP884252DOA-B
<b>Abbreviation</b>	PCMP-E102
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FIH2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	VTKSPNS IPELVKHIDS DLIRNAHKVF DEIPELDVIS ATAVIGRFVK ESRHVEASQA FKRLCLGIR PNEFTFGTVI GSSTTSRDVK LGKQLHCYAL KMGLASNVFV GSAVLNCYVK LSTLTDARRC FDDTRDPNVV SITNLISGYL KKHEFEEALS LFRAMPERSV VTWNAVIGGF SQTGRNEEAV NTFVDMLREG VVIPNESTFP CAITAINIA SHGAGKSIHA CAIKFLGKRF NVFVWNSLIS FYSKCGNMED SLLAFNKLEE EQRNIVSWNS MIWGYAHNGR GEEAVAMFEK MVKDTNLRPN NVTILGVLFA CNHAGLIQEG YMYFNKAVND YDDPNLLELE HYACMVDMLS RSGRFKEAEE LIKSMPLDPG IGFWKALLGG CQIHSNKRLA KLAASKILEL DPRDVSSYVM LSNAYSAMEN WQNVSLIRRK MKETGLKRFT GCSWIEVRDQ IRVFNADKN NELKDEVYRM LALVSQHLEE NECWKDL
<b>Source</b>	E.coli
<b>Target Names</b>	PCMP-E102
<b>Protein Names</b>	Recommended name: Pentatricopeptide repeat-containing protein At5g42450, mitochondrial
<b>Expression Region</b>	24-517
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.