



Recombinant Arabidopsis thaliana Putative F-box/FBD/LRR-repeat protein At5g56810 (At5g56810)

Product Code	CSB-MP884265DOA
Abbreviation	At5g56810
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FJT2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEAARSGIED VTSPDRISQL PNDLLFRILS LIPVSDAMST SLLSKRWKSV WKMLPTLVYN ENSCSNIGSL GFDQFCGRSL QLHEAPLLKT LTLELRKQTD SLDSSIFPNI HSTLLEFSIK STGYPVYYST ISFPNLDVF QTLVVLLKLQG NICLDVVDSP VCFQSLKSLY LTCVNFENEE SFSKLLSACP VLEDLFLQRL CSVGRFLFSI SVPSLQRLTY TKEQAYYSND EAILEITAPS LKHLNIFDRV GVFSFIEDMP KLVEASVRVK LSKNEKLPKV LTSVEHLSLD LYPSMVFHLD DRFISKQLLH LKLDIYDNFQ SNLLLSLLKD LPNLQSLKLN HSHPSYNVED QPCSVSEPS VPECLSFHLE TFQWIGYAGT FEEIAAAVYV LKNARCLKNA TISLYSRGTE NGLMMIKELE SMSKGSIMCQ LLVKF
Source	Mammalian cell
Target Names	At5g56810
Protein Names	Recommended name: Putative F-box/FBD/LRR-repeat protein At5g56810
Expression Region	1-435
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.