



Recombinant Arabidopsis thaliana Pantoate--beta-alanine ligase (At5g48840)

Product Code	CSB-MP884271DOA
Abbreviation	At5g48840
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FKB3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEPEVIRDKD SMRKWSRAMR SQGKTIGLVP TMGYLHEGHL SLVRQSLALT DVTVVSIVVN PGQFSPTEDL STYPSDFSGD LTKLAALSGG KVVVFNPKNL YDYGGETKKI NDGGGNGGRV VSCVEEGGLG HETWIRVERL EKGFCGKSRP VFFRGVATIV TKLFNIVEPD VALFGKKDYQ QWRRIQRMVR DLNFGIEIVG SDIAREKDGL AMSSRNVRLS DEERQRALSI SRSLAMAKAS VAEGKTNCAE LKDMIIQQVV GSAGRVDYVE IVDQETLEGV EEIKSGVVIC VAAWFGTVRL IDNIEINVSL
Source	Mammalian cell
Target Names	At5g48840
Protein Names	Recommended name: Pantoate--beta-alanine ligase EC= 6.3.2.1 Alternative name(s): Pantoate-activating enzyme Pantothenate synthetase
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.