



Recombinant Arabidopsis thaliana EID1-like F-box protein 2 (EDL2)

Product Code	CSB-MP884293DOA
Abbreviation	EDL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FLZ8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MIIAKQYRCI HSATCHCTKG HLSEEVFLFM VQHLNWNPNV IATLSCVCKW FDDLAKRLLW KEFCRARAPK MMSDLQSSGS HSVDGSWRAL GKLLIYCSGS SKGGLFNDVQ ISGHFVHRTR FSRTSGRSFL PPQCRTDDIL YVSDPCEHLD QGEGDGLGFF RGIFKSFSMS KVRKLLIKKG TPFHPTEVCP YCKAKLWSML QAKMIPQSAS CRLGAYEDSI EYYVCLNGHM LGVCTLLPLS DSEGASEFQ
Source	Mammalian cell
Target Names	EDL2
Protein Names	Recommended name: EID1-like F-box protein 2
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.