



# Recombinant Arabidopsis thaliana Polyamine oxidase 1 (PAO1)

<b>Product Code</b>	CSB-MP884307DOA
<b>Abbreviation</b>	PAO1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FNA2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSTASVIIIG AGISGISA AK VLVENGVEDV LILEATDRIG GRIHKQNF GD VPVELGAGWI AGVGGKESNP VWELASRFNL RTCFSDYTNA RFNIYDRSGK IFPTGIASDS YKKA VDSAIL KLKSLEAQCS GQVAEEAPSS PKTPIELAID FILHDFEMAE VEPISTYVDF GEREFLVADE RGYECLLYKM AEEFLVTSHG NILDYRLKLN QVVREVQQSR NGVVVKTEDG SVYEANYVIV SASIGVLQSD LLSFQPLLPR WKTEAIQKCD VMVYTKIFLK FPQCFWPCGP GQEFFIYAHE QRGYFTFWQH MENAYPGSNI LVVTLTNEQS KRVEAQSDQE TMKEAMSVLR DMFGATIPYA TDILVPRWWN NRFQRGSYSN YPMISDNQLL QNIKAPVGRI FFTGEHTSEK FSGYVHGGYL AGIDTSKSL L EEMKQSLLLQ PLLAFTESLT LTHQKPNNSQ IYTNVKFISG TS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PAO1
<b>Protein Names</b>	Recommended name: Polyamine oxidase 1 Short name= AtPAO1 EC= 1.5.3.16 Alternative name(s): N(1)-acetylpolyamine oxidase Spermine oxidase
<b>Expression Region</b>	1-472
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.