



# Recombinant Arabidopsis thaliana Zinc finger CCCH domain-containing protein 63 (ZFWD2)

<b>Product Code</b>	CSB-BP884314DOA
<b>Abbreviation</b>	ZFWD2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FNZ1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDFDLNGGNK RVFNRLGGGG GSTRPMAPTD TRQKVCFHWR AGRCNRSPCP YLHRELPGPG PGQQQGGPGYT NKRVAEESGF AGPSHRRGPG FNGNSSSSSWG RFGGNRTVTK TEKVCNFWVD GNCTYGDKCR YLHCWSKGES FALLTQLDGH EKLVSIALP SGSDKLYTGS KDETLRVWDC ASGQCTGVK LGGEIGCVLS EGPWLLVGMP NLVKAWNIET NADQSLSGPV GQVYSLVGT DLLFAGTQDG SILAWRYNAA TNCFEPSASL TGHTLAVVTL YVGANRLYSG SMDKTIKVWS LDNLQCIQTL TDHSSVMSL ICWDQFLLSC SLDNTVKIWA AIEGGNLEVT YTHKEEHGVL ALCGVHDAEA KPVLLCACND NTLRLYDLPS LGLFIRFTER GKIFAKQEIR AIQIGPGGIF FTGDGTGQVK VWKWCTEPTA ALP
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZFWD2
<b>Protein Names</b>	Recommended name: Zinc finger CCCH domain-containing protein 63 Short name= AtC3H63 Alternative name(s): Zinc finger CCCH domain and WD40 repeat-containing protein 2
<b>Expression Region</b>	1-443
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.