



Recombinant *Oryza sativa* subsp. *japonica* REF/SRPP-like protein Os05g0151300/LOC_Os05g05940 (Os05g0151300, LOC_Os05g05940)

Product Code	CSB-MP884320OFG
Abbreviation	Os05g0151300, LOC_Os05g05940
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FRA7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MADSGSDAPI SNRPEEEVTV EKTPEMEAAA EEERLRYLEF VQQAAAQVLV LAAAAYAYAK QGAGPLRPGV DHVEGTVKAV VGPVYDRFHG VPLDLLKFLD RKVGESVQEL DRRVPPVVKE APGLARSAAA EVRQAGLVGT ATGLAKSAIA RAEPRARDLY TRYEPVAERK AAEAWAALNR LPLVPSVTRA VLPAAASLSA RYNTAVADGA KRGSATVATYL PLVPTERLSR VFGYPLADAA TSPAPEMQPI PSQ
Source	Mammalian cell
Target Names	Os05g0151300
Protein Names	Recommended name: REF/SRPP-like protein Os05g0151300/LOC_Os05g05940
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.