



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 2 (Os01g0192000, LOC_Os01g09620)

Product Code	CSB-MP884327OFG
Abbreviation	Os01g0192000, LOC_Os01g09620
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FU27
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	<pre> MMMMGEGVSS VPPWSHLPVS GVDVLGGGGG GGDEMTPYVI AALRDYLPAN DVGVGADEEE EAAAMAAVD AYACDEFMY EFKVRRRCARG RSHDWTECPF AHPGEKARRR DPRKYHYSGT ACPDFRKGGC KRGDACEYAH GVFECLWHPA RYRTQPCKDG TACRRRVCF AHTPDQLRVL PAQQSSPRSV ASSPLAESYD GSPLRRQAFE SYLTKTIMSS SPTSTLMSPP KSPPESEPL SPDGAAAIRR GSWPGVGSPV NDVLASFRQL RLNKVKSSPS GGWSYPSSSA VYGSPKAATG LYSLPTTPLA STATVTTASS FMPNLEPLDL GLIGDEEPVQ RVESGRALRE KVFERLSRDG AISGDATAFA TAGVGLDVDW VSDLIN </pre>
Source	Mammalian cell
Target Names	Os01g0192000
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 2 Short name= OsC3H2 Alternative name(s): Protein DELAY OF THE ONSET OF SENESCENCE Short name= OsDOS
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.