



Recombinant Arabidopsis thaliana 12-oxophytodienoate reductase 3 (OPR3)

Product Code	CSB-MP884329DOA
Abbreviation	OPR3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FUP0
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MTAAQGNSNE TLFSSYKMGR FDLSHRVVLA PMTRCRLNG VPNAALAEYY AQRTPGGFL ISEGMTVSPG SAGFPHVPGI YSDEQVEAWK QVVEAVHAKG GFIFCQLWHV GRASHAVYQP NGGSPISSTN KPISENRWRV LLPDGSHVKY PKPRALEASE IPRVVEDYCL SALNAIRAGF DGIEIHGAHG YLIDQFLKDG INDRTDQYGG SIANRCRFLK QVVEGVVSAI GASKVGVVRS PAIDHLDATD SDPLSLGLAV VGMLNKLQGV NGSKLAYLHV TQPRYHAYGQ TESGRQGSDE EEAKLMKSLR MAYNGTFMSS GGFNKELGMQ AVQQGDADLV SYGRLFIANP DLVSRFKIDG ELNKYNRKTf YTQDPVVGyT DYPFLAPFSR L
Source	Mammalian cell
Target Names	OPR3
Protein Names	Recommended name: 12-oxophytodienoate reductase 3 EC= 1.3.1.42 Alternative name(s): 12-oxophytodienoate-10,11-reductase 3 Short name= AtOPR3 Short name= OPDA-reductase 3 Delayed dehiscence 1
Expression Region	1-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.