



# Recombinant Arabidopsis thaliana Shaggy-related protein kinase delta (ASK4)

<b>Product Code</b>	CSB-MP884335DOA
<b>Abbreviation</b>	ASK4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9FVS6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESHLGNGVG SSRSKNTKN TSSVDWLSR DMLEMKIRDK TEADEERDSE PDIIDGVGAE PGHVITTTLP GRNGQSRQTV SYIAEHVVG TGSFGMV FQAK CRETGEVVAI KVLQDKRYK NRELQIMQML DHPNVVCLKH SFYSRTENEE VYLNLVLEFV PETVNRTARS YSRMNQLMPL IYVKLYTYQI CRGLAYLHNC CGLCHRDIP QNLLVNPHTH QLKICDFGSA KVLVKGEPNI SYICSRYYRA PELIFGATEY TTAIDIWSTG CVMAELLGQ PLFPGESGVD QLVEIHKVLG TPTREEIKCM NPNYTEFKFP QIKPHPWHKV FQKRLPPEAV DLLCRFFQYS PNLRCTAVEA CIHPFFDEL RDPNARLPNGR PLPPLFNFKP QELSGIPPET VDRLVPEHAR KQNHFMALHS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ASK4
<b>Protein Names</b>	Recommended name: Shaggy-related protein kinase delta EC= 2.7.11.1 Alternative name(s): ASK-delta
<b>Expression Region</b>	1-420
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.