



Recombinant Arabidopsis thaliana F-box only protein 6 (FBX6)

Product Code	CSB-BP884359DOA
Abbreviation	FBX6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9FZK1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEEELAMLRQ LIGQLQELLH NGSPPPSSSS SSLSSSSPSF LVLHHPQYQN GWCLPCIEDT SADDCCDIVM AGGKRPGIFK MLETVKPPVK RTRKERTQGK SCTEVDEISG NMDQEIWQEF PQDLFEDVVS RLPMATFFQF RAVCRKWNAL IDSDSFRCF TELPQTIPWF YTITHENVNS GQVYDPSLKK WHHPPIPALP KKSIVLPMAS AGGLVCFLDI GHRNFYVSNP LTKSFRELP RSKVWSRVA VGMTLNGNST SHGYKVLWVG CEGEYEVYDS LSNVWTKRGT IPSNIKLPVL LNFKSQPVAI HSTLYFMLTD PEGILSYDMV SGKWKQFIIP GPPDLSDHTL AACGERLMLV GLLTKNAATC VCIWELQKMT LLWKEVDRMP NIWCLEFYGK HIRMNCLGNK GCLILLSLRS RQMNRLITYN AVTREWTKVP GCTVPRGRKR LWIACGTAFH PSPTARA
Source	Baculovirus
Target Names	FBX6
Protein Names	Recommended name: F-box only protein 6
Expression Region	1-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.