



Recombinant Enterobacteria phage AR1 Receptor-recognizing protein (38)

Product Code	CSB-MP884362EEAN
Abbreviation	38
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9G0B4
Product Type	Recombinant Protein
Immunogen Species	Escherichia phage AR1 (Bacteriophage AR1)
Purity	>85% (SDS-PAGE)
Sequence	MAVTGPWVGS SAVVNTGQNW MVGAAQRLRM GAPFWMSNMI GRSVEVIHTL GADHNFNGQW FRDRCFEAGS APIVFNITGD LVSYSRDVPL FFMYGDTPNE YVQLNIHGVT MYGRGGNGWA ACAIGASDGG VCIQNDIGGR LRINNGGAIA GGGGGGGGYS QANNWAGKYV CGGGGGRPFG LGGNNGARWP GGNASLTSPG AGGNTGTRY Y AGGGGEVGQP GQYANPGASY STPPTSPGAA VAGSAPTWQN VGAIYGPRV
Source	Mammalian cell
Target Names	38
Protein Names	Recommended name: Receptor-recognizing protein Alternative name(s): Protein Gp38
Expression Region	1-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.