




# Recombinant Human Protein Wnt-10a (WNT10A)

<b>Product Code</b>	CSB-YP884424HU
<b>Relevance</b>	Ligand for members of the frizzled family of seven transmembrane receptors. Probable developmental protein. May be a signaling molecule important in CNS development. Is likely to signal over only few cell diameters.
<b>Abbreviation</b>	Recombinant Human WNT10A protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9GZT5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MPRSAPNDILDLRLPPEPVLNANTVCLTLPGLSRRQMEVCVRHPDVAASAIQGI QIAIHECQHQRDQRWNCSSLETRNKIPYESPIFSRGFRESAFAYAIAAAGVVH AVSNACALGKLKACGCDASRRGDEEAFRRKLHRLQLDALQRGKGLSHGVPE HPALPTASPLQDSWEWGGCSPDMGFGERFSKDFLDSREPHRDIHARMRLH NNRVGRQAVMENMRRKCKCHGTSGSCQLKTCWQVTPEFRTVGALLRSRFH RATLIRPHNRNGGQLEPGPAGAPSPAPGAPGPRRRASPADLVYFEKSPDFCE REPR LDSAGTVGRLCNKSSAGSDGCGSMCCGRGHNILRQTRSERCHCRFH WCCFVVCEEWRITEWVSVCK
<b>Research Area</b>	Cancer
<b>Source</b>	Yeast
<b>Target Names</b>	WNT10A
<b>Protein Names</b>	Recommended name: Protein Wnt-10a
<b>Expression Region</b>	36-417aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	44.4kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.