



Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 2 (CRISPLD2)

Product Code	CSB-EP884430HU-B
Abbreviation	CRISPLD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H0B8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	YLLPNVTL LEELLSKYQH NESHSRVRRRA IPREDKKEIL MLHNKLRGQV QPQASNMEYM TWDDELEKSA AAWASQCIWE HGPTSLLVSI GQNLGAHWGR YRSPGFHVQS WYDEVKDYTY PYPSECNPWC PERCSGPMCT HYTQIVWATT NKIGCAVNTC RKMTVWGEVW ENAVYFVCNY SPKGNWIGEA PYKNGRPCSE CPPSYGGSCR NNLCYREETY TPKPETDEM EVETAPIPEE NHVWLQPRVM RPTKPKKTSV VNYMTQVVRD DTKMKDRCKG STCNRYQCPA GCLNHKAKIF GTLFYESSSS ICRAAIHYGI LDDKGGGLVDI TRNGKVPFFV KSERHGVQSL SKYKPSSSFV VSKVKVQDLQ CYTTVAQLCP FEKPATHCPR IHCPAHCKDE PSYWAPVFGT NIYADTSSIC KTAVHAGVIS NESGGDQDVM PVDKKTQYVG SLRNGVQSES LGTPRDGKAF RIFAVRQ
Source	E.coli
Target Names	CRISPLD2
Protein Names	Recommended name: Cysteine-rich secretory protein LCCL domain-containing 2 Alternative name(s): Cysteine-rich secretory protein 11 Short name= CRISP-11 LCCL domain-containing cysteine-rich secretory protein 2
Expression Region	23-497
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.