



Recombinant Human Testis-specific Y-encoded-like protein 1 (TSPYL1)

Product Code	CSB-YP884436HU
Abbreviation	TSPYL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H0U9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSGLDGVKRT TPLQTHSIII SDQVPSDQDA HQYLRLRDQS EATQVMAEPG EGGSETVALP PPPPSEEGGV PQDAAGRGGT PQIRVVGGRG HVAIKAGQEE GQPPAEGGLAA ASVVMAADRS LKKGVQGGEK ALEICGAQRS ASELTAGAEA EAEEVKTGKC ATVSAAVAER ESAEVVKEGL AEKEVMEEQM EVEEQPPEGE EIEVAEEDRL EEEAREEEEGP WPLHEALRMD PLEAIQLELD TVNAQADRAF QQLEHKFGRM RRHYLERRNY IIQNIPGFWM TAFRNHPQLS AMIRGQDAEM LRYITNLEVK ELRHPRTGCK FKFFFRRNPY FRNKLIVKEY EVRSSGRVVS LSTPIIWRRG HEPQSFIRRN QDLICFFFTW FSDHSLPESD KIAEIIKEDL WPNPLQYYLL REGVRRARRR PLREPVEIPR PFGFQSG
Source	Yeast
Target Names	TSPYL1
Protein Names	Recommended name: Testis-specific Y-encoded-like protein 1 Short name= TSPY-like protein 1
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is found in the nucleolus and is similar to that of a family of genes on the Y-chromosome. This gene is intronless. Defects in this gene are a cause of sudden infant death with dysgenesis of the testes syndrome (SIDDT).
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.