



Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 1 (CRISPLD1)

Product Code	CSB-EP884457HU-B
Abbreviation	CRISPLD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H336
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVVPNAT LLEKLEKYM DEDGEWWIAK QRGKRAITDN DMQSILDLHN KLRSQVYPTA SNMEYMTWDV ELERSAESWA ESCLWEHGPA SLLPSIGQNL GAHWGRYRPP TFHVQSWYDE VKDFSYPYEH ECNPYCPFRC SGPVCTHYTQ VVWATSNRIG CAINLCHNMN IWGQIWPKAV YLVCNYSKPG NWWGHAPYKH GRPCSACPPS FGGGCRENLK YKEGSDRYYP PREEETNEIE RQQSQVHDTH VRTRSDDSSR NEVISAQQMS QIVSCEVRLR DQCKGTTNCR YECPAGCLDS KAKVIGSVHY EMQSSICRAA IHYGIIDNDG GWVDITRQGR KHYFIKSNRN GIQTIGKYQS ANSFTVSKVT VQAVTCETTV EQLCPFHKPA SHCPRVYCPR NCMQANPHYA RVIGTRVYSD LSSICRAAVH AGVVRNHGGY VDVMPVDKRRK TYIASFQNGI FSESLQNPPG GKAFRVFAVV
Source	E.coli
Target Names	CRISPLD1
Protein Names	Recommended name: Cysteine-rich secretory protein LCCL domain-containing 1 Alternative name(s): CocoaCrisp Cysteine-rich secretory protein 10 Short name= CRISP-10 LCCL domain-containing cysteine-rich secretory protein 1 Trypsin
Expression Region	24-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.