



Recombinant Human Cysteine-rich secretory protein LCCL domain-containing 1 (CRISPLD1)

| | |
|--------------------------|---|
| Product Code | CSB-MP884457HU |
| Abbreviation | CRISPLD1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9H336 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MVVPNAT LLEKLEKYM DEDGEWWIAK QRGKRAITDN DMQSILDLHN KLRSQVYPTA SNMEYMTWDV ELERSAESWA ESCLWEHGPA SLLPSIGQNL GAHWGRYRPP TFHVQSWYDE VKDFSYPYEH ECNPYCPFRC SGPVCTHYTQ VVWATSNRIG CAINLCHNMN IWGQIWPKAV YLVCNYSPPKG NWWGHAPYKH GRPCSACPPS FGGGCRENLC YKEGSDRYYP PREEETNEIE RQQSQVHDTH VRTRSDDSSR NEVISAQQMS QIVSCEVRLR DQCKGTTCSR YECPAGCLDS KAKVIGSVHY EMQSSICRAA IHYGIIDNDG GWVDITRQGR KHYFIKSNRN GIQTIGKYQS ANSFTVSKVT VQAVTCETTV EQLCPFHKPA SHCPRVYCPR NCMQANPHYA RVIGTRVYSD LSSICRAAVH AGVVRNHGGY VDVMPVDKRR TYIASFQNGI FSESLQNPPG GKAFRVFAVV |
| Source | Mammalian cell |
| Target Names | CRISPLD1 |
| Protein Names | Recommended name: Cysteine-rich secretory protein LCCL domain-containing 1 Alternative name(s): CocoaCrisp Cysteine-rich secretory protein 10 Short name= CRISP-10 LCCL domain-containing cysteine-rich secretory protein 1 Trypsin |
| Expression Region | 24-500 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.