



Recombinant Human Probable tRNA threonylcarbamoyladenosine biosynthesis protein OSGEPL1 (OSGEPL1)

Product Code	CSB-BP884466HU
Abbreviation	OSGEPL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H4B0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	P GTLFLHKIVL GIETSCDDTA AAVVDETGNV LGEAIHSQTE VHLKTGGIVP PAAQQLHREN IQRIVQEALS ASGVSPSDL S AIATTIKPGL ALSLGVGLSF SLQLVGQLKK PFIPIHHMEA HALTIRLTNK VEFPLVLLI SGGHCLLALV QGVSDFLLLG KSLDIAPGDM LDKVARRLSL IKHPECSTMS GGKAIEHLAK QGNRFHFDIK PPLHHAKNCD FSFTGLQHVT DKIIIMKKEKE EGIEKGQILS SAADIAATVQ HTMACHLVKR THRILFCKQ RDLLPQNNAV LVASGGVASN FYIRRALEIL TNATQCTLLC PPPRLCTDNG IMIAWNGIER LRAGLGILHD IEGIRYEPKC PLGVDISKEV GEASIKVPQL KMEI
Source	Baculovirus
Target Names	OSGEPL1
Protein Names	Recommended name: Probable tRNA threonylcarbamoyladenosine biosynthesis protein OSGEPL1 Alternative name(s): t(6)A37 threonylcarbamoyladenosine biosynthesis protein OSGEPL1
Expression Region	30-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.