



# Recombinant Human Probable tRNA threonylcarbamoyladenosine biosynthesis protein OSGEPL1 (OSGEPL1)

<b>Product Code</b>	CSB-MP884466HU
<b>Abbreviation</b>	OSGEPL1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H4B0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	P GTLFLHKIVL GIETSCDDTA AAVVDETGNV LGEAIHSQTE VHLKTGGIVP PAAQQLHREN IQRIVQEALS ASGVSPSDL S AIATTIKPGL ALSLGVGLSF SLQLVGQLKK PFIPIHHMEA HALTIRLTNK VEFPLVLLI SGGHCLLALV QGVSDFLLLG KSLDIAPGDM LDKVARRLSL IKHPECSTMS GGKAIEHLAK QGNRFHFDIK PPLHHAKNCD FSFTGLQHVT DKIIMKKEKE EGIEKGQILS SAADIAATVQ HTMACHLVKR THRILFCKQ RDLLPQNNAV LVASGGVASN FYIRRALEIL TNATQCTLLC PPPRLCTDNG IMIAWNGIER LRAGLGILHD IEGIRYEPKC PLGVDISKEV GEASIKVPQL KMEI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	OSGEPL1
<b>Protein Names</b>	Recommended name: Probable tRNA threonylcarbamoyladenosine biosynthesis protein OSGEPL1 Alternative name(s): t(6)A37 threonylcarbamoyladenosine biosynthesis protein OSGEPL1
<b>Expression Region</b>	30-414
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.