



Recombinant Human Differentially expressed in FDCP 6 homolog (DEF6)

Product Code	CSB-MP884467HU
Abbreviation	DEF6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H4E7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MALRKELLKS IWYAFTALDV EKSGKVSQSQ LKVLSHNLYT VLHIPHDPVA LEEHRDDDD GPVSSQGYMP YLNKYILDKV EEGAFVKEHF DELCWTLTAK KNYRADSNGN SMLSNQDAFR LWCLFNFLSE DKYPLIMVPD EVEYLLKKVL SSMSLEVSLG ELEELLAQEA QVAQTTGGLS VWQFLELFNS GRCLRGVGRD TLSMAIHEVY QELIQDVLKQ GYLWKRGLR RNWAERWFQL QPSCLCYFGS EECKEKGRII PLDAHCCVEV LPDRDGKRCM FCVKTANRTY EMSASDTRQR QEWTAAIQMA IRLQAEGKTS LHKDLKQKRR EQREQRERRR AAKEEELLRL QQLQEEKERK LQELELLQEA QRQAERLLQE EEERRRSQHR ELQQALEGQL REAEQARASM QAEMELKEEE AARQRRIKE LEEMQQRLQE ALQLEVKARR DEESVRIAQT RLLEEEEEKL KQLMQLKEEQ ERYIERAQQE KEELQQEMAQ QSRSLQQAQQ QLEEVQRNRQ RADEDVEAAQ RKLQRQASTNV KHWNVQMNRL MHPIEPGDKR PVTSSSFSGF QPPLLAHRDS SLKRLTRWGS QGNRTPSPNS NEQQKSLNGG DEAPAPASTP QEDKLDPAPE N
Source	Mammalian cell
Target Names	DEF6
Protein Names	Recommended name: Differentially expressed in FDCP 6 homolog Short name= DEF-6 Alternative name(s): IRF4-binding protein
Expression Region	1-631
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.