



Recombinant Human RIB43A-like with coiled-coils protein 2 (RIBC2)

Product Code	CSB-EP884469HU
Abbreviation	RIBC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H4K1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRQNDKIMCI LENRKKRDRK NLCRAINDFQ QSFQKPETRR EFDLSDPLAL KKDLPARQSD NDVRNTISGM QKFMGEDLNF HERKKFQEEQ NREWSLQQQR EWKNARAEQK CAEALYTETR LQFDETAKHL QKLESTTRKA VCASVKDFNK SQAIESVERK KQEKKQEVED NLAEITNLLR GDLLSENPPQ AASSFGPHRV VPDRWKGMTQ EQLEQIRLVQ KQQIQEKLRL QEEKRQRDL WDRRRIQGAR ATLLFERQQW RRQRDLRRAL DSSNLSLAKE QHLQKKYMNE VYTNQPTGDY FTQFNTGSR
Source	E.coli
Target Names	RIBC2
Protein Names	Recommended name: RIB43A-like with coiled-coils protein 2
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.