



Recombinant Human EH domain-containing protein 1 (EHD1)

Product Code	CSB-BP884470HU
Abbreviation	EHD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H4M9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFSWVSKDAR RKKEPELFQT VAEGLRQLYA QKLLPLEEHY RFHEFHSPAL EDADFDNKPM VLLVGQYSTG KTTFIRHLIE QDFPGMRIGP EPTTDSFIAV MHGPTEGVVP GNALVVDPRR PFRKLNAFGN AFLNRFMCAQ LPNPVLDSIS IIDTPGILSG EKQRISRGYD FAAVLEWFAE RVDRILLFD AHKLDISDEF SEVIKALKNH EDKIRVVLNK ADQIETQQLM RVYGALMWSL GKIINTPEVV RVYIGSFWSH PLLIPDNRKL FEAEEQDLFK DIQSLPRNAA LRKLNLIKR ARLAKVHAYI ISSLKKEMPN VFGKESKKKE LVNNLGEIYQ KIEREHQISP GDFPSLRKMQ ELLQTQDFSK FQALKPKLLD TVDDMLANDI ARLMVMVRQE ESLMPSQVVK GGAFDGTMNG PFGHGYGEGA GEGIDDVEWV VGKDKPTYDE IFYTLSPVNG KITGANAKKE MVKSKLPNTV LGKIWKLADV DKDGLLDDEE FALANHLIKV KLEGHELPAD LPPHLVPPSK RRHE
Source	Baculovirus
Target Names	EHD1
Protein Names	Recommended name: EH domain-containing protein 1 Alternative name(s): PAST homolog 1 Short name= hPAST1 Testilin
Expression Region	1-534
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene belongs to a highly conserved gene family encoding EPS15 homology (EH) domain-containing proteins. The protein-binding EH domain was first noted in EPS15, a substrate for the epidermal growth factor receptor. The EH domain has been shown to be an important motif in proteins involved in protein-protein interactions and in intracellular sorting. This protein is thought to play a role in the endocytosis of IGF1 receptors.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.