



Recombinant Human mRNA (2'-O-methyladenosine-N (6)-)-methyltransferase (PCIF1)

Product Code	CSB-BP884471HU
Abbreviation	PCIF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H4Z3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MANENHGSPR EEASLLSHSP GTSNQSQPCS PKPIRLVQDL PEELVHAGWE KCWSRRENRP YYFNRFNTQS LWEMPVLGQH DVISDPLGLN ATPLPQDSSL VETPPAENKP RKRQLSEEQP SGNGVKKPKI EIPVPTTGQS VPSSPSIPGT PTLKMWGTSP EDKQQAALLR PTEVYWDLDI QTNAVIKHRG PSEVLPPHPE VELLSQLIL KLRQHYRELC QQREGIEPPR ESFNRWMLER KVVDKGSDPL LPSNCEPVVS PSMFREIMND IPIRLSRIKF REEAKRLLFK YAEAARRLIE SRSASPDSRK VVKWNVEDTF SWLRKDHSAS KEDYMDRLEH LRRQCGPHVS AAKDSVEGI CSKIYHISLE YVKRIREKHL AILKENNISE EVEAPEVEPR LVYCYPVRLA VSAPPMPSTE MHMENNVCV RYKGMVKVS RNYFSKLWLL YRYSCIDDSA FERFLPRVWC LLRRYQMMFG VGLYEGTGLQ GSLPVHVFEA LHRLFGVSFE CFASPLNCYF RQYCSAFPDT DGYFGSRGPC LDFAPLSGSF EANPPFCEEL MDAMVSHFER LLESSPEPLS FIVFIPEWRE PPTPALTRME QSRFKRHQLI LPAFEHEYRS GSQHICKKEE MHYKAVHNTA VLFLQNDPGF AKWAPTPERL QELSAAYRQS GRSHSSGSSS SSSSEAKDRD SGREQGPSRE PHPT
Source	Baculovirus
Target Names	PCIF1
Protein Names	Recommended name: Phosphorylated CTD-interacting factor 1
Expression Region	1-704
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.