



# Recombinant Human mRNA (2'-O-methyladenosine-N (6)-)-methyltransferase (PCIF1)

<b>Product Code</b>	CSB-EP884471HU-B
<b>Abbreviation</b>	PCIF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H4Z3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MANENHGSPR EEASLLSHSP GTSNQSQPCS PKPIRLVQDL PEELVHAGWE KCWSRRENRP YYFNRFNTQS LWEMPVLGQH DVISDPLGLN ATPLPQDSSL VETPPAENKP RKRQLSEEQP SGNGVKKPKI EIPVTPTGQS VPSSPSIPGT PTLKMWGTSP EDKQQAALLR PTEVYWDLDI QTNAVIKHRG PSEVLPPHPE VELLRSQLIL KLRQHYRELC QQREGIEPPR ESFNRWMLER KVDKGS DPL LPSNCEPVVS PSMFREIMND IPIRLSRIKF REEAKRLLFK YAEAARRLIE SRSASPD SRK VVKWNVEDTF SWLRKDHSAS KEDYMDRLEH LRRQCGPHVS AAKDSVEGI CSKIYHISLE YVKRIREKHL AILKENNISE EVEAPEVEPR LVYCYPVRLA VSAPPMP SVE MHMENNVC I RYKGEMVKVS RNYFSKLWLL YRYSCIDDSA FERFLPRVWC LLRRYQMMFG VGLYEGTGLQ GSLPVHVFEA LHRLFGVSFE CFASPLNCYF RQYCSAFPDT DGYFGSRGPC LDFAPLSGSF EANPPFCEEL MDAMVSHFER LLESSPEPLS FIVFIPEWRE PPTPALTRME QSRFKRHQLI LPAFEHEYRS GSQHICKKEE MHYKAVHNTA VLFLQNDPGF AKWAPTPERL QELSAAYRQS GRSHSSGSSS SSSSEAKDRD SGREQGPSRE PHPT
<b>Source</b>	E.coli
<b>Target Names</b>	PCIF1
<b>Protein Names</b>	Recommended name: Phosphorylated CTD-interacting factor 1
<b>Expression Region</b>	1-704
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.