



Recombinant Human Putative uncharacterized protein FLJ22184

Product Code	CSB-MP884478HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6K5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPPLETQSSL APPSLQTPPA SLTTPPLENL PSLAPPPLQT ASAPLTTPHL ETPPCPAPCP LQAPPSPLTT PPPETPSSIA TPPPQAPPAL ASPPLQGLPS PPLSPLATPP PQAPPALALP PLQAPPSPPA SPPLSPLATP SPQAPNALAV HLLQAPFSPS PAPPVQAPFS PPASPPVSPS ATPPSQAPPS LAAPPLQVPP SPPASPPMSP SATPPPQAPP PLAAPPLQVP PSPPASPPMS PSATPPPRVP PLLAAPPLQV PPSPASLPM SPLAKPPPQA PPALATPPLQ ALPSPPASFP GQAPFSPSAS LPMSPLATPP PQAPPVLAAP LLQVPPSPPA SPTLQAPRRP PTPGPDTSVS GPRLTLALAP GPPPPPSRSP SSTLSGPDLA GHSSSATSTP EELRGYDSGP EGGAAASPPP DAELAACHPA AWSRGPAPPP AFRGAPGAPL PWPPATGPGS ADGLCTIYET EGPESATPAP GALDPGPS PG TSGGKAAAGA GAGASSRSPK QARLGELPLG ALQASVVQHL LSRTLLAAA EGAAGGSGGG PGGAEGGGVT GGARAALSDA ELGRWAELLS PLDESRASIT SVTSFSPDDV ASPQGDWTVV EVETFH
Source	Mammalian cell
Target Names	PRR36
Protein Names	Recommended name: Putative uncharacterized protein FLJ22184
Expression Region	1-616
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.