



Recombinant Human Uncharacterized protein C19orf44 (C19orf44)

Product Code	CSB-EP884480HU-B
Abbreviation	C19orf44
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H6X5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MASARKASRP MRDVFGDFSD VSLEDSTMEE IRNFQISRNL TKIAPGHSRF LKRNQTLDK HLLKKNPVL GSGPRLASCR PPTASRIRA NAALMKLAQL ETRIMNRKLQ RNLSDTESDS MTADAGLPKR ADRILSGGAL ELASQNTDKT SQNQARELPV TENNAQNAKV SRFLKKKQAP VENISPEAPA GKERTLQTPK QKEPARTFDS PDSDEEEMKV LLSLMDSSR EKNTNQGFSS ANVSEEEERK LFSVPSQLRA FTVPSVELSS AKPSQTSHP TSLAADRTLH STRSRADYPQ SHVSSDTASH TPSVSITGAF SNSVSLKMGH VKLVSSPGRS EAETVDEPVS EGADDSLDEF RINILSLDGL APAVSENSDL EQEEESAQRQ KTAGKIFRAE ASTGQDAPRQ AQARSWASQG KAASAEGDES EVSEHLSASS ASAIQQDSTS SMQPPSEAPM VNTVSSAYSE DFENSPSLTA SEPTAHSKES LDRTLDALSE SSSSVKTDLP QTAESRKKSG RHVTRVLVKD TAVQTPDPAF TYEWTKVASM AAMGPALGGA YVDPTPIANH VISADAIEAL TAYSPAVLAL HDVLKQQLSL TQQFIQASRH LHASLLRSLD ADSFHYHTLE EAKEYIRCHR PAPLTMEDAL EEVNKEL
Source	E.coli
Target Names	C19orf44
Protein Names	Recommended name: Uncharacterized protein C19orf44
Expression Region	1-657
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.