



# Recombinant Human Uncharacterized protein C19orf44 (C19orf44)

<b>Product Code</b>	CSB-MP884480HU
<b>Abbreviation</b>	C19orf44
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H6X5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASARKASRP MRDVFGDFSD VSLEDSTMEE IRNFQISRNL TKIAPGHSRF LKRNQTLDK HLLKKNPVL GSGPRLASCR PPTASRIRA NAALMKLAQL ETRIMNRKLQ RNLSDTESDS MTADAGLPKR ADRILSGGAL ELASQNTDKT SQNQARELPV TENNAQNAKV SRFLKKKQAP VENISPEAPA GKERTLQTPK QKEPARTFDS PDSDEEEMKV LLSLMDSSR EKNTNQGFSS ANVSEEEERK LFSVPSQLRA FTVPSVELSS AKPSQTSHP TSLAADRTLH STRSRADYPQ SHVSSDTASH TPSVSITGAF SNSVSLKMGH VKLVSSPGRS EAETVDEPVS EGADDSLDEF RINILSLDGL APAVSENSDL EQEEESAQRQ KTAGKIFRAE ASTGQDAPRQ AQARSWASQG KAASAEGDES EVSEHLSASS ASAIQQDSTS SMQPPSEAPM VNTVSSAYSE DFENSPSLTA SEPTAHSKES LDRTLDALSE SSSSVKTDLP QTAESRKKSG RHVTRVLVKD TAVQTPDPAF TYEWTKVASM AAMGPALGGA YVDPTPIANH VISADAIEAL TAYSPAVLAL HDVLKQQLSL TQQFIQASRH LHASLLRSLD ADSFHYHTLE EAKEYIRCHR PAPLTMEDAL EEVNKEL
<b>Source</b>	Mammalian cell
<b>Target Names</b>	C19orf44
<b>Protein Names</b>	Recommended name: Uncharacterized protein C19orf44
<b>Expression Region</b>	1-657
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.