



Recombinant Human Uncharacterized protein KIAA0226-like (KIAA0226L)

Product Code	CSB-MP884481HU
Abbreviation	KIAA0226L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H714
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVSQSTVRQD SPVEPWEGIS DHSGIIDGSP RLLNTDHPPC QLDIRLMRHK AVWINPQDVQ QQPQDLQSQV PAAGNSGTHF VTDAASPSGP SPSCCLGDSLA ETTLSEDTTD SVGSASPHGS SEKSSSFSLs STEVHMVRPG YSHRVSLPTS PGILATSPYP ETDSAFFEPS HL TSAADEGA VQVSRRTISS NSFSPEVFVL PVDVEKENAH FYVADMIISA MEKMKCNILS QQQTESWSKE VSGLLGSDQP DSEMTFDTNI KQESGSSTSS YSGYEGCAVL QVSPVTETRT YHDVKEICKC DVDEFVILEL GDFNDITETC SCSCSSSKSV TYEPDFNSAE LLAKELYRVF QKCWILSVVN SQLAGSLSAA GSIVVNEECV RKDFESSMNV VQEIKFKSRI RGTEDWAPPR FQIIFNIHPP LKRDLVVAQ NFFCAGCGTP VEPKFVKRLR YCEYLGKYFC DCCHSYAESC IPARILMMWD FKKYYVSNFS KQLLDSIWHQ PIFNLLSIGQ SLYAKAKELD RVKEIQEQLF HIKKLLKTCR FANSALKEFE QVPGHLTDEL HLFSLDLVR IKKGLLAPLL KDILKASLAH VAGCELCQGK GFICEFCQNT TVIFPFQTAT CRRCSACRAC FHKQCFQSSE CPRCARITAR RKLLESVASA AT
Source	Mammalian cell
Target Names	RUBCNL
Protein Names	Recommended name: Uncharacterized protein KIAA0226-like
Expression Region	1-662
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.