



Recombinant Human YY1-associated protein 1 (YY1AP1)

Product Code	CSB-EP884488HU-B
Abbreviation	YY1AP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H869
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEEASRSAA ATNPGSRLTR WPPPKREGS AVDPGKRRSL AATPSSSLPC TLIALGLRHE KEANELMEDL FETFQDEMGE SNMEDDGPEE EERVAEPQAN FNTPQALRFE ELLANLLNEQ HQIAKELFEQ LKMKKPSAKQ QKEVEKVKPQ CKEVHQTLLI DPAQRKRLQQ QMQQHVVLLT QIHLLATCNP NLNPEASSTR ICKELGTFA QSSIALHHQY NPKFQTLFQP CNLMGAMQLI EDFSTHVSID CSPHKTVKKT ANEFPCLPKQ VAWILATSKV FMYPELLPVC SLKAKNPQDK ILFTKAEDNK YLLTCKTARQ LTVRIKLNLM NRAPDNIIFK YKKTQQLPVL GKCCEEIQPH QWKPPIEREE HRLPFWLKAS LPSIQEELRH MADGAREVGN MTGTTEINSD QGLEKDENSEL GSETRYPLLL PKGVVLKLP VADRFPPKAW RQKRSSVLKP LLIQPSPSLQ PSFNPGKTPA QSTHSEAPPS KMLRIPHI QPATVLQTVG VPPPLGVSGG ESFESPAALP AMPPEARTSF PLSESQTLLS SAPVPKVMMP SPASSMFRKP YVRRRPSKRR GARAFRCIKP APVIHPASVI FTVPATTVKI VSLGGGCNMI QPVNAAVAQS PQTIPIATLL VNPTSFPCL NQPLVASSVS PLIVSGNSVN LPIPSTPEDK AHMNVDIACA VADGENAFQG LEPKLEPQEL SPLSATVFPK VEHSPGPPPV DKQCQEGLSE NSAYRWTVVK TEEGRQALEP LPQGIQESLN NSSPGDLEEV VKMEPEDATE EISGFL
Source	E.coli
Target Names	YY1AP1
Protein Names	Recommended name: YY1-associated protein 1 Alternative name(s): Hepatocellular carcinoma susceptibility protein Hepatocellular carcinoma-associated protein 2
Expression Region	1-796
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The encoded gene product presumably interacts with YY1 protein; however, its



exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.