



Recombinant Human Cyclin N-terminal domain-containing protein 2 (CNTD2)

Product Code	CSB-MP884491HU
Abbreviation	CNTD2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H8S5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLVRGRDQGS GSRLGPIVRR WAPRPSPLQS LAASLDAEPS SAAVPDGFFPA GPTVSPRRLA RPPGLEEALS ALGLQGEREY AGDIFAEVMV CRVLPLRALP RAVTPEMRAL VVDWLWQVHE YLGLAGDTLY LAVHLLDSYL SAGRVRLHRL QLLGVACLFV ACKMEECVLP EPAFLCLLSA DSFSRAELLR AERRILSRDL FRLHHPGPLL CLGLLAALAG SSPQVMLLAT YFLELSLLEA EAAGWEPGRR AAAALSLAHR LLDGAGSRLQ PELYRCSLGG GSVWGHRSFR DLPSWSFLRS RRMRDNY
Source	Mammalian cell
Target Names	CNTD2
Protein Names	Recommended name: Cyclin N-terminal domain-containing protein 2
Expression Region	1-307
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.