



Recombinant Human Pleckstrin homology domain-containing family F member 2 (PLEKHF2)

Product Code	CSB-MP884492HU
Abbreviation	PLEKHF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H8W4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVDRLANSEA NTRRISIVEN CFGAAGQPLT IPGRVLIGEG VLTKLCRKKP KARQFFLFND ILVYGNIVIQ KKKYNKQHII PLENTIDSI KDEGDLRNGW LIKTPTKSFA VYAATATEKS EWMNHINKCV TDLLSKSGKT PSNEHAAVWV PDSEATVCMR CQKAKFTPVN RRHHCRCXGF VVCGPCSEKR FLLPSQSSKP VRICDFCYDL LSAGDMATCQ PARSDSYSQS LKSPLNDMSD DDDDDSSD
Source	Mammalian cell
Target Names	PLEKHF2
Protein Names	Recommended name: Pleckstrin homology domain-containing family F member 2 Short name= PH domain-containing family F member 2 Alternative name(s): PH and FYVE domain-containing protein 2 Phafin-2 Zinc finger FYVE domain-containing pr
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.