



# Recombinant Human Ankyrin repeat and zinc finger domain-containing protein 1 (ANKZF1)

<b>Product Code</b>	CSB-MP884493HU
<b>Abbreviation</b>	ANKZF1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H8Y5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSPAPDAAPA PASISLFDLS ADAPVFQGLS LVSHAPGEAL ARAPRTSCSG SGERESPERK LLQGPMDisE KLCSTCDQT FQNHQEQREH YKLDWHRFNL KQRLKDKPLL SALDFEKQSS TGDLSISGS ESDSASEED LQTLDRERAT FEKLSRPPGF YPHRVLFQNA QGQFLYAYRC VLGPHQDPPE EAELLQNLQ SRGPRDCVVL MAAAGHFAGA IFQGREVVTH KTFHRYTVRA KRGTAAQLRD ARGGPPSHSAG ANLRRYNEAT LYKDVRDLLA GPSWAKALEE AGTILLRAPR SGRSLFFGGK GAPLQRGDPR LWDIPLATRR PTFQELQRVL HKLTTLHVYE EDPREAVRLH SPQTHWKTVR EERKKPTEEE IRKICRDEKE ALGQNEESPK QGSSEGEDG FQVELELVEL TVGTLDLCES EVLPKRRRRK RNKKEKSRDQ EAGAHRTLLQ QTQEEEPSTQ SSQAVAAPLG PLLDEAKAPG QPELWNALLA ACRAGDVGVL KLQLAPSPAD PRVLSLLSAP LGSGGFTLLH AAAAAGRGSV VRLLEAGAD PTVQDSRARP PYTVAADKST RNEFRRFMEK NPDAYDYNKA QVPGPLTPEM EARQATRKRE QKAARRQREE QQQRQQEQEE REREEQRRFA ALSDREKRAL AAERRLAAQL GAPTSPIPDS AIVNTRRCWS CGASLQGLTP FHYLDFSFCS TRCLQDHRRQ AGRPSS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ANKZF1
<b>Protein Names</b>	Recommended name: Ankyrin repeat and zinc finger domain-containing protein 1 Alternative name(s): Zinc finger protein 744
<b>Expression Region</b>	1-726
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.