



Recombinant Human Ankyrin repeat and zinc finger domain-containing protein 1 (ANKZF1)

Product Code	CSB-YP884493HU
Abbreviation	ANKZF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H8Y5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSPAPDAAPA PASISLFDLS ADAPVFQGLS LVSHAPGEAL ARAPRTSCSG SGERESPERK LLQGPMDisE KLCSTCDQT FQNHQEQREH YKLDWHRFNL KQRLKDKPLL SALDFEKQSS TGDLSISGS ESDSASEED LQTLDRERAT FEKLSRPPGF YPHRVLFQNA QGQFLYAYRC VLGPHQDPPE EAELLQNLQ SRGPRDCVVL MAAAGHFAGA IFQGREVVTH KTFHRYTVRA KRGTAAQLRD ARGGPHSAG ANLRRYNEAT LYKDVRDLLA GPSWAKALEE AGTILLRAPR SGRSLFFGGK GAPLQRGDPR LWDIPLATRR PTFQELQRVL HKLTTLHVYE EDPREAVRLH SPQTHWKTVR EERKKPTEEE IRKICRDEKE ALGQNEESPK QGSSEGEDG FQVELELVEL TVGTLDLCES EVLPKRRRRK RNKKEKSRDQ EAGAHRTLLQ QTQEEEPSTQ SSQAVAAPLG PLLDEAKAPG QPELWNALLA ACRAGDVGVL KLQLAPSPAD PRVLSLLSAP LGSGGFTLLH AAAAAGRGSV VRLLEAGAD PTVQDSRARP PYTVAADKST RNEFRFMMEK NPDAYDYNKA QVPGPLTPEM EARQATRKRE QKAARRQREE QQQRQQEQEE REREEQRRFA ALSDREKRAL AAERRLAAQL GAPTSPIPDS AIVNTRRCWS CGASLQGLTP FHYLDFSFCS TRCLQDHRRQ AGRPSS
Source	Yeast
Target Names	ANKZF1
Protein Names	Recommended name: Ankyrin repeat and zinc finger domain-containing protein 1 Alternative name(s): Zinc finger protein 744
Expression Region	1-726
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.