



Recombinant Human CCR4-NOT transcription complex subunit 10 (CNOT10)

Product Code	CSB-MP884495HU
Abbreviation	CNOT10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H9A5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AADKPADQG AEKHEGTGQS SGITDQEKEL STNAFQAFTS GNYDACLQHL ACLQDINKDD YKIILNTAVA EFFKSNQTTT DNLRQTLNQL KNQVHSAVEE MDGLDDVENS MLYYNQAVIL YHLRQYTEAI SVGEKLYQFI EPFEEKFAQA VCFLLDVLYI LTYQAEKALH LLAVLEKMIS QGNNNKNGKN ETGNNNNKDG SNHKAESGAL IEAAKSKIHQ YKVRAYIQMK SLKACKREIK SVMNTAGNSA PSLFLKSNFE YLRGNYRKAV KLLNSSNIAE HPGFMKTGEC LRCMFWNNLG CIHFAMSKHN LGIFYFKKAL QENDNVCAQL SAGSTDPGKK FSGRPMCTLL TNKRYELLYN CGIQLLHIGR PLAAFECLIE AVQVYHANPR LWLRLAECCI AANKGTSEQE TKGLPSKKGI VQSIVGQGYH RKIVLASQSI QNTVYNDGQS SAIPVASMEF AAICLRNALL LLPEEQQDPK QENGAKNSNQ LGGNTESSES SETCSSKSHD GDKFIPAPPS SPLRKQELLEN LKCSILACSA YVALALGDNL MALNHADKLL QQPKLSGSLK FLGHLYAAEA LISLDRISDA ITHLNPENVT DVSLGISSNE QDQGSDDKGEN EAMESSGKRA PQCYPSSVNS ARTVMLFNLG SAYCLRSEYD KARKCLHQAA SMIHPKEVPP EAILLAVYLE LQNGNTQLAL QIIKRNQLLP AVKTHSEVRK KPVFQPVHPI QPIQMPAFTT VQRK
Source	Mammalian cell
Target Names	CNOT10
Protein Names	Recommended name: CCR4-NOT transcription complex subunit 10
Expression Region	2-744
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.