



# Recombinant Human CCR4-NOT transcription complex subunit 10 (CNOT10)

<b>Product Code</b>	CSB-YP884495HU
<b>Abbreviation</b>	CNOT10
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9H9A5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	AADKPADQG AEKHEGTGQS SGITDQEKEL STNAFQAFTS GNYDACLQHL ACLQDINKDD YKIILNTAVA EFFKSNQTTT DNLRQTLNQL KNQVHSAVEE MDGLDDVENS MLYYNQAVIL YHLRQYTEAI SVGEKLYQFI EPFEEKFAQA VCFLLVDLYI LTYQAEKALH LLAVLEKMIS QGNNNKNGKN ETGNNNNKDG SNHKAESGAL IEAAKSKIHQ YKVRAYIQMK SLKACKREIK SVMNTAGNSA PSLFLKSNFE YLRGNYRKAV KLLNSSNIAE HPGFMKTGEC LRCMFWNNLG CIHFAMSKHN LGIFYFKKAL QENDNVCAQL SAGSTDPGKK FSGRPMCTLL TNKRYELLYN CGIQLLHIGR PLAAFECLIE AVQVYHANPR LWLRLAECCI AANKGTSEQE TKGLPSKKG I VQSIVGQGYH RKIVLASQSI QNTVYNDGQS SAIPVASM EF AAICLRNALL LLPEEQQDPK QENGA KNSNQ LGGNTESSES SETCSSKSHD GDKFIPAPPS SPLRKQEL EN LKCSILACSA YVALALGDNL MALNHADKLL QQP KLSGSLK FLGHL YAAEA LISLDRISDA ITHLNPENVT DVSLGISSNE QDQGS DKGEN EAMESSGKRA PQCYPSSVNS ARTVMLFNLG SAYCLRSEYD KARKCLHQAA SMIHPKEVPP EAILLAVYLE LQNGNTQLAL QIIKR NQLLP AVKTHSEVRK KPVFQPVHPI QPIQMPAFTT VQRK
<b>Source</b>	Yeast
<b>Target Names</b>	CNOT10
<b>Protein Names</b>	Recommended name: CCR4-NOT transcription complex subunit 10
<b>Expression Region</b>	2-744
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.